SHORT REPORT Open Access

Suitability of individual and bulk milk samples to investigate the humoral immune response to lumpy skin disease vaccination by ELISA



Milovan Milovanović¹, Vesna Milićević², Sonja Radojičić¹, Miroslav Valčić¹, Bernd Hoffmann³ and Klaas Dietze^{3*}

Abstract

Background: The detection of antibodies against capripoxvirus has become easier with a commercially available ELISA validated for serum and plasma. In order to explore its suitability for immunological investigations on alternative samples, this study targeted milk as sample matrix available through non-invasive sampling.

Methods: Samples for this study were collected from dairy cows vaccinated against LSD in an area without reported LSD virus circulation. Paired serum and milk (individual and bulk) samples were tested by ELISA without and with modifications of the sample incubation time for the milk samples. For the evaluation of the test specificity, 352 milk samples from a milk repository in Germany were used as negative control. Receiver operating characteristic analysis was performed for determination of the Youden index and determination of the most suitable cut-off value for maximum specificity.

Results: From 154 analyzed serum samples from Serbia, 75 were detected as positive in the ELISA. Sensitivity and specificity of the ELISA test for milk samples reached values of 88 to 91% using Youden criteria. A cut-off of 10 was determined aiming for maximum specificity. This cut-off value was used for further analysis. Using the protocol for serum, out of 154 milk samples, 38 were detected as positive, number of positive detected milk samples increase up to 48 with modified protocol. Milk samples from Germany reacted negative, except two samples that had borderline results using modified protocol. Significant statistical difference (p < 0.05) was observed between two incubation protocols. The detection of LSD-specific antibodies from bulk milk samples (pools of 2–10 individuals) came along with a reduced sensitivity over the sample of individual animals.

Conclusions: Results show that the detection of capripoxvirus specific antibodies in milk samples using the commercially available ELISA from IDvet is feasible and can represent a helpful tool for LSDV monitoring programs.

Keywords: Lumpy skin disease, ELISA, Milk, Herd screening, Surveillance, Serology, Non-invasive

Introduction

Lumpy skin disease (LSD) is a vector-borne disease of cattle caused by lumpy skin disease virus (LSDV) which belongs to genus *Capripoxvirus* [1, 2]. Clinically, the disease can manifest in a wide spectrum ranging from severe acute to sub-acute and in-apparent forms. Typical signs observed are fever, enlargement of lymph nodes,

nasal discharge, and firm skin nodules [3]. The control of LSD is mainly based on mass vaccination of the susceptible cattle population with live attenuated capripoxvirus vaccine [4, 5].

Recent studies suggest that LSDV vaccination stimulates equally cell mediated and humoral immunity [6–8]. The humoral immune response is of paramount practical importance to obtain information on the immune status of animals post infection or vaccination. The detection of antibodies against LSDV is possible starting one to 2 weeks post vaccination, increases gradually

Full list of author information is available at the end of the article



^{*} Correspondence: klaas.dietze@fli.de

³Friedrich-Loeffler-Institut, Südufer 10, D-17493 Greifswald-Insel Riems, Germany

until 35 days to 12 weeks post vaccination, and is described to last until 40 weeks post vaccination [6–8]. Standard serological methods like virus neutralization test (VNT), agar gel immunodiffusion, indirect fluorescence antibody test (IFAT) and Western blot are described [9, 10].

These tests are expensive and time consuming therefore limiting their use for fast serological screening of cattle populations. The only test validated to date by the OIE is the VNT, coming along with high specificity and good sensitivity, but reduced high-throughput application. As performing VNT includes the handling of live capripoxvirus, its application can face additional restrictions [11]. The IFAT comes along with the disadvantage of increased cross reactivity with bovine papular stomatitis virus and other poxvirus antibodies.

For efficient immunological investigations at large scale, Enzyme-linked Immunosorbent Assay (ELISA) has been found to be more suitable compared to the abovementioned serological methods. Few studies on the development of ELISA tests for LSD-specific antibody detection by using recombinant P32 [12], recombinant two virion core protein of sheeppox virus [13], and inactivated sheeppox virus as coating antigen for antibody detection [14] have been published. A double antigen ELISA from the IDvet® is currently the only commercially available kit for detection of capripox specific antibodies applied in field studies for seromonitoring [15]. All of these ELISA have been developed for plasma or serum as sample matrix to be used.

Apart from blood, the detection of specific antibodies is also possible from other matrices such as milk. Cow milk contains three major classes of immunoglobulin's (Ig): IgG, IgM and IgA [16]. Dominant class of immunoglobulins in milk, comprising about 65% of the total, is IgG. The concentration of IgG varies between serum and mammary secretions. Highest concentrations can be found in colostrum (32-212 mg/ml), followed by serum (25.0 mg/ml) and finally in milk with about 0.72 mg/ml [17]. As for the colostrum, a steady gradual decrease of antibody levels is described each hour [18]. The concentration of antibodies in regular milk is much lower than that in colostrum and depends on factors like clinical or subclinical mastitis, age, breed, feeding system and the stage of lactation. In addition, primiparous cows are described to have colostrum and milk with overall lower levels of immunoglobulins compared to multiparous cows [18, 19].

Using milk for the detection of antibodies by ELISA has been found suitable targeting antibodies against viruses such as bovine viral diarrhea virus (BVDV) [20], bluetongue virus [21, 22] and bovine alpha herpesvirus 1 (BHV-1) [23]. As milk samples are non-invasive and cheaper to collect compared to serum and plasma, this

method has been widely used in mass screening activities targeting either individual animals or bulk milk samples for the determination of the immunological status at animal or herd level respectively.

In this study, we investigated the suitability of milk, as individual animal and bulk sample to detect LSD-specific antibodies using a commercially available capripox ELISA designed for serum and plasma samples.

Material and methods

Sample material

For the purpose of this study, serum and milk samples were collected from 154 lactating dairy cows from farms in the Kraljevo municipality in Serbia. This area had no evidence of previous LSD field virus circulation. All cows included in this study were vaccinated in 2016 with a Neethling vaccine strain (Onderstepoort Biological Products, South Africa) and re-vaccinated in 2017 with BOVIVAX LSD-N Neethling vaccine (M.C.I. Sante Animale, Morocco). During sampling, all animals were clinically examined. Animals showing any generic clinical signs of disease were excluded from study. Blood samples were collected using blood colleting tubes without anticoagulant BD Vacutainertm4055269 (Belliver Industrial Estate, UK) via coccygeal venipuncture. Extracting serum from blood was done by centrifugation at 2000 RPM for 20 min. Milk samples were collected in 15 ml sterile plastic tubes, individual milk samples by hand milking from all teats and bulk milk from milk tank of the farm. Milk samples were conserved using 1% sodium-azide. All obtained samples were aliquoted in 1.5 ml centrifugal tubes and stored at -20 °C until further examination. Repository samples of cattle from Germany (n = 352) were included in this study as negative control panel since Germany is LSD-free and the vaccination against LSDV is forbidden.

Elisa

For the detection of LSD-specific antibodies in serum and milk samples a commercial ELISA test for serum and plasma (ID Screen* Capripox Double Antigen Multi-species ELISA, IDvet, Montpellier, France) was used. Investigation of serum and milk (individual and bulk) samples was performed according to manufacture protocol for serum and plasma. Additionally, investigation of milk (individual and bulk) samples was performed with modifications of sample incubation time, from 90 min at $+21\,^{\circ}\text{C}$ to overnight incubation at $+4\,^{\circ}\text{C}$. All milk samples were centrifuged at 10000 rpm for 10 min; fat supernatant was removed and only the liquid fraction was used for analysis.

Milovanović et al. Virology Journal (2020) 17:28 Page 3 of 7

Statistical analysis

Statistical analysis was performed using RStudio® Version 1.1.456 (Boston, USA) with the packages Optimal Cutpoints Version 1.1–3 and reshape2 1.4.3 from the CRAN repository. Receiver operating characteristic (ROC) analysis was performed for determination of the Youden index and determination of the most suitable cut-off value for maximum specificity for milk samples using serum ELISA results with as reference. The McNemar test was performed to determine the statistical significance between results of the two ELISA protocols used for the milk samples.

Results

Sampling and clinical examination

No clinical mastitis or any generic clinical signs of disease were observed during sample collection. During the 2017 sampling period, no new LSD outbreak was reported in Serbia. In total 154 individual and 38 bulk milk samples were collected. The 38 obtained bulk milk samples comprised of milk from two to ten individuals depending on the overall farm size and the number of animals in lactation.

Antibody detection

From 154 investigated serum samples, 75 were detected as positive with ELISA test. The results obtained from the serum served as reference for the ROC analysis evaluating the two different protocols used on the milk samples to determine the respective cut-off values for the criteria Youden and maximum specificity. The results are summarized in Table 1 and Fig. 1.

Following these cut-off determinations, all further analysis of the milk samples used the cut-off of ≥ 10 as determined through the maximum specificity criterion.

From the tested 154 milk samples, samples from 38 animals tested positive when applying the proposed manufacturer protocol for serum and plasma with a 90 min incubation time. The number of positive samples from the same sample set increased to 48 when using the modified overnight incubation without creating false positive results in reference to the serum. The relationship of the results on milk samples using two different protocols, one treating milk samples identical to the serum samples and the other using an overnight incubation time is depicted in Fig. 2. The results between the two incubation protocols are significantly different (p < 0.05), with overnight incubation obtaining higher S/P% values.

The overall agreement between investigated serum and milk samples with both protocols is shown in Table 2, stratified additionally by the level of positivity in the serum ELISA. Higher levels of agreement are achieved the higher the serum ELISA values are.

The sensitivity of ELISA for the bulk milk sample is however reduced over the testing of the individual animal whilst the specificity remained high. No false positive bulk milk samples appeared in relation to obtained results of individual milk samples composing it. Bulk prevalence had a close relation with bulk milk results showing better performance of bulk milk samples, which had higher prevalence. The composition and results of 38 bulk milk samples are summarized in Fig. 3, for details of the bulk milk composition refer to Additional file 1.

Except two samples from Germany, which had borderline result using overnight incubation, the rest of the samples reacted negative.


Discussion

In the light of the LSD expansion in recent years, tools to determine the immune status of animals and populations are on demand to facilitate control and eradication

Table 1 Comparative ROC analysis of milk samples with 90 min and overnight incubation time using Youden and maximum specificity as criteria

Criterion Incubation	Youden		Maximum specificity	
	90 min	overnight	90 min	overnight
Area under ROC curve (ACU)	0.961 (0.942, 0.98)	0.932 (0.887, 0.977)	0.959 (0.939, 0.978)	0.92 (0.87, 0.971)
cut-off	1.308	3.086	9.425	9.785
Se	0.895	0.882	0.513	0.632
Sp	0.891	0.914	1.0	1.0
PPV	0.591	0.644	1.0	1.0
NPV	0.980	0.978	0.921	0.939
DLR. Positive	8.186	10.222	Inf	Inf
DLR. Negative	0.118	0.13	0.487	0.368
FP	47.0	37.0	0.0	0.0
FN	8.0	9.0	37.0	28.0
Optimal criterion	0.785	0.795		

Milovanović et al. Virology Journal (2020) 17:28 Page 4 of 7



efforts. The obtained results from this study provide the proof of principle, that milk can be used as sample for the detection of antibodies against LSDV both at the individual animal and herd level.

Applying overnight incubation for milk samples at +4 °C in ELISA tests has been described as suitable

for the improved detection of specific antibodies for different viruses [20, 23, 24]. In this study, where no further modifications on neither the milk sample nor the ELISA kit have been applied, the increased incubation time confirms to be suitable to achieve higher test sensitivity.

Milovanović et al. Virology Journal (2020) 17:28 Page 5 of 7




Fig. 2 Detected S/P% values of serum samples and their respective milk samples following two incubation protocols. Cut-off for serum at 30 as provided by manufacturer and for milk set at 10 according to ROC analysis criterion maximum specificity

The determination of the optimal cut-offs for will depend on the actual purpose of the testing. In this study higher cut-off value aiming for maximum specificity was taken for further analysis, accommodating a situation where false positive samples would have a higher negative impact than false negative, as all animals were exposed to LSDV through vaccination. For a final confirmation, serum samples could be collected and analyzed by ELISA or other serological tests recommended by OIE [1] if needed.

In general, observed mismatches of detected positive serum and negative milk samples together with lower S/P% values in milk samples, could be attributed to lower concentration of immunoglobulins in milk which is multifactorial [18, 19]. Improvement of the ELISA test performance for the detection of LSD-specific antibodies

in milk samples can probably be addressed by purification of immunoglobulins using commercially available kits or by precipitation of proteins using ammonium sulfate [25]. On the other side, Klintevall et al. [24] showed in their study that cows that suffered of clinical mastitis may give false positive ELISA results when investigating milk. Cows included in this study did not suffer from clinical mastitis and therefore this aspect was not observed.

The results obtained in this study were not put in relation to the exact date of vaccination of the animals. It is well described however that the levels of detectable antibodies in the serum vary and will increase between one and 12 weeks post vaccination, after which period antibody titer will slowly start to decrease [6]. Similar developments of titers are likely to be seen in the milk

Table 2 Comparison of detected positive and negative milk samples to detected S/P% values of positive serum samples

Positive serum samples		Individual milk samples				
		90 min incubation		Overnight incubation		
S/P% value	Number	Positive	Negative	Positive	Negative	
≥201	28	25 (89.29%)	3 (10.71%)	27 (96.43%)	1 (3.57%)	
101–200	24	9 (37.50%)	15 (62.50%)	16 (66.67%)	8 (33.33%)	
30–100	23	4 (17.39%)	19 (82.61%)	5 (21.74%)	18 (78.26%)	

Milovanović et al. Virology Journal (2020) 17:28 Page 6 of 7




Fig. 3 Bulk milk results stratified over the within bulk prevalence determined through individual milk testing; overnight incubation and cut-off value ≥10

especially after vaccination campaigns stop. In case of a LSDV re-emergence, the detection of positive milk samples will imply that animals were in contact with the virus provoking antibody synthesis.

One additional step forward in exploring the diagnostic advantages of milk samples was the detection of LSD-specific antibodies in pooled milk samples representing milk bulks typical for Serbia. Benefits of using bulk milk samples have been described previously for fast preliminary immunological investigation at herd level in response of determination of seroprevalence against BVDV and BHV-1 [26]. The detection of antibodies from bulk milk mainly depends of the number of positive animals included in bulk milk and the concentration of immunoglobulins in positive milk included in the bulk milk. Even though the number of animals per bulk in this study was not too big, detection of positive bulk milk samples was possible starting with bulk prevalence of 25%. As expected, higher bulk milk prevalence delivered better results by detecting all positive bulk milk samples which had prevalence higher that 50%. According to obtained bulk milk results and described LSD morbidity rate of 45% [27], determination of LSDV seroprevalence using bulk milk samples can be suitable. Multiple testing of bulk milk can be additionally recommended on a wider scale for screening purposes which will provide robust, but fast herd level results. In addition, these results should be validated for settings with larger animal numbers contributing to the bulk milk.

Conclusion

The commercially available ELISA kit from IDvet for the detection of capripoxvirus specific antibodies is in principle suitable to be used on milk samples, from

individual animals as well as pooled milk samples of small bulks. Cut-off values will need to be specified according to the purpose of testing. Additional modifications either on the ELISA kit or through antibody enrichment steps can lead to improve the sensitivity and specificity and will make this approach a time and cost efficient, non-invasive monitoring for LSD occurrence or LSD vaccination in the field.

Supplementary information

Supplementary information accompanies this paper at https://doi.org/10.1186/s12985-020-01298-x.

Additional file 1. Table providing an overview on bulk milk sample composition and the respective ELISA results.

Abbreviations

BHV-1: Bovine alpha herpesvirus 1; BVDV: Bovine viral diarrhea virus; ELISA: Enzyme-linked Immunosorbent Assay; IFAT: Immunofluorescent antibody test; Ig: Immunoglobulin; LSD: Lumpy skin disease; LSDV: Lumpy skin disease virus; npv: Predictive negative value; ppv: Predictive positive value; ROC: Receiver operating characteristic; S/P %: Percentage of positivity compare to the positive control; Se: Sensitivity; Sp: Specificity; VNT: Virus neutralization test

Acknowledgements

Many thanks to German Academics Exchange Service (DAAD) section ST21 – East Central Europe, East South Europe, and Turkey. The authors thank the farmers in Serbia for their participation as well as Christian Korthase, Cindy Krenz, Kim Lea Molle and Karin Pinger for the excellent technical assistance in the laboratory.

Authors' contributions

MM, VM, SR, BH and KD planned and designed the study. MM collected the sample material. MM and VM processed samples after collection and prepared them for analysis. MM carried out the sample analysis. MM, KD and BH interpreted results and drafted the manuscript. SR, VM and MV critically revised the manuscript. All authors have read the manuscript and approved submission of the manuscript.

Funding

Collection and shipment of sample material were supported by the Ministry of Education, Science and Technological Development of the Republic of Serbia (project grant numbers TR31088 and TR37015).

Availability of data and materials

The datasets used and/or analysed during the current study are available from the corresponding author and/or the first author on reasonable request.

Ethics approval and consent to participate

The study was approved by the Ministry of Agriculture, Forest and Water Management, Veterinary Directorate, Serbia. The respective reference number is 323/07/5369/2018/05/2.

Consent for publication

Not applicable.

Competing interests

The authors declare that they have no competing interests.

Author details

¹Department of Infectious Diseases of Animals and Diseases of Bees, University of Belgrade, Faculty of Veterinary Medicine, Blvd. Oslobodjenja 18, Belgrade 11000, Serbia. ²Virology Department, Institute of Veterinary Medicine of Serbia, Vojvode Toze 14, Belgrade 11000, Serbia. ³Friedrich-Loeffler-Institut, Südufer 10, D-17493 Greifswald-Insel Riems, Germany

Received: 12 August 2019 Accepted: 18 February 2020 Published online: 05 March 2020

References

- OIE. Lumpy Skin Disease technical disease card 2017 [Available from: www.oie.int
- 2. Tuppurainen ES, Oura CA. Review: lumpy skin disease: an emerging threat to Europe, the Middle East and Asia. Transbound Emerg Dis. 2012;59(1):40–8.
- Babiuk S, Bowden TR, Parkyn G, Dalman B, Manning L, Neufeld J, et al. Quantification of lumpy skin disease virus following experimental infection in cattle. Transbound Emerg Dis. 2008;55(7):299–307.
- Tuppurainen ES, Pearson CR, Bachanek-Bankowska K, Knowles NJ, Amareen S, Frost L, et al. Characterization of sheep pox virus vaccine for cattle against lumpy skin disease virus. Antivir Res. 2014;109:1–6.
- Gari G, Abie G, Gizaw D, Wubete A, Kidane M, Asgedom H, et al. Evaluation of the safety, immunogenicity and efficacy of three capripoxvirus vaccine strains against lumpy skin disease virus. Vaccine. 2015;33(28):3256–61.
- Abdelwahab MG, Khafagy HA, Moustafa AM, Saad MA. Evaluation of Humoral and cell-mediated immunity of lumpy skin disease vaccine prepared from local strainin calves and its related to maternal immunity. J Am Sci. 2016;12(10):8.
- Norian R, Afzal Ahangaran N, Hr V, Azadmehr A. Evaluation of Humoral and cell-mediated immunity of two Capripoxvirus vaccine strains against lumpy skin disease virus. Iranian J Virol. 2016;10(4):1–11.
- Varshovi HR, Norian R, Azadmehr A, Afzal Ahangaran N. Immune response characteristics of Capri pox virus vaccines following emergency vaccination of cattle against lumpy skin disease virus. Iranian J Vet Sci Technol. 2018;9(2):7.
- OIE Manual of Diagnostic Tests and Vaccines for Terrestrial Animals 2018. Chapter 3.4.12. Paris 2018. (https://www.oie.int).
- Gari G, Biteau-Coroller F, LeGoff C, Caufour P, Roger F. Evaluation of indirect fluorescent antibody test (IFAT) for the diagnosis and screening of lumpy skin disease using Bayesian method. Vet Microbiol. 2008;129(3–4):269–80.
- Tuppurainen ESM, Venter EH, Shisler JL, Gari G, Mekonnen GA, Juleff N, et al. Review: Capripoxvirus diseases: current status and opportunities for control. Transbound Emerg Dis. 2017;64(3):729–45.
- Carn VM, Kitching RP, Hammond JM, Chand P. Use of a recombinant antigen in an indirect ELISA for detecting bovine antibody to capripoxvirus. J Virol Methods. 1994;49(3):285–94.
- 13. Bowden TR, Coupar BE, Babiuk SL, White JR, Boyd V, Duch CJ, et al. Detection of antibodies specific for sheeppox and goatpox viruses using

- recombinant capripoxvirus antigens in an indirect enzyme-linked immunosorbent assay. J Virol Methods. 2009;161(1):19–29.
- Babiuk S, Wallace DB, Smith SJ, Bowden TR, Dalman B, Parkyn G, et al. Detection of antibodies against Capripoxviruses using an inactivated Sheeppox virus ELISA. Transbound Emerg Dis. 2009;56(4):132–41.
- Milovanović M, Dietze K, Milićević V, Radojičić S, Valčić M, Moritz T, et al. Humoral immune response to repeated lumpy skin disease virus vaccination and performance of serological tests. BMC Vet Res. 2019:15(1):80.
- Korhonen H, Marnila P, Gill HS. Milk immunoglobulins and complement factors. Br J Nutr. 2000;84(Suppl 1):575–80.
- 17. Larson BL. Immunoglobulins of the mammary secretions. New York: Elsevier;
- Morin DE, Nelson SV, Reid ED, Nagy DW, Dahl GE, Constable PD. Effect of colostral volume, interval between calving and first milking, and photoperiod on colostral IgG concentrations in dairy cows. J Am Vet Med Assoc. 2010;237(4):420–8.
- Krol J, Litwinczuk Z, Brodziak A, Barlowska J. Lactoferrin, lysozyme and immunoglobulin G content in milk of four breeds of cows managed under intensive production system. Pol J Vet Sci. 2010;13(2):357–61.
- Beaudeau F, Belloc C, Seegers H, Assie S, Sellal E, Joly A. Evaluation of a blocking ELISA for the detection of bovine viral diarrhoea virus (BVDV) antibodies in serum and milk. Vet Microbiol. 2001;80(4):329–37.
- Mars MH, van Maanen C, Vellema P, Kramps JA, van Rijn PA. Evaluation of an indirect ELISA for detection of antibodies in bulk milk against bluetongue virus infections in the Netherlands. Vet Microbiol. 2010;146(3): 209–14.
- Kramps JA, van Maanen K, Mars MH, Popma JK, van Rijn PA. Validation of a commercial ELISA for the detection of bluetongue virus (BTV)-specific antibodies in individual milk samples of Dutch dairy cows. Vet Microbiol. 2008;130(1):80–7.
- Nylin B, Stroger U, Ronsholt L. A retrospective evaluation of a bovine Herpesvirus-1 (BHV-1) antibody ELISA on bulk-tank milk samples for classification of the BHV-1 status of Danish dairy herds. Prev Vet Med. 1999;47(1–2):91–105.
- Klintevall K, Näslund K, Svedlund G, Hajdu L, Linde N, Klingeborn B.
 Evaluation of an indirect ELISA for the detection of antibodies to bovine leukaemia virus in milk and serum. J Virol Methods. 1991;33(3):319–33.
- Wingfield PT. Protein Precipitation Using Ammonium Sulfate. Curr Protoc Protein Sci. 2016;84(1):A.3F.1–9.
- Stahl K, Rivera H, Vagsholm I, Moreno-Lopez J. Bulk milk testing for antibody seroprevalences to BVDV and BHV-1 in a rural region of Peru. Prev Vet Med. 2002;56(3):193–202.
- 27. Coetzer JAW, Tuppurainen ES, Babiuk S, Wallace DB. Lumpy skin disease. 2018 [Available from: www.anipedia.org.

Publisher's Note

Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

Ready to submit your research? Choose BMC and benefit from:

- fast, convenient online submission
- thorough peer review by experienced researchers in your field
- rapid publication on acceptance
- support for research data, including large and complex data types
- gold Open Access which fosters wider collaboration and increased citations
- maximum visibility for your research: over 100M website views per year

At BMC, research is always in progress.

Learn more biomedcentral.com/submissions

